ALSATOM SU - MPC

ELETTROBISTURI PER CHIRURGIA MONOPOLARE E BIPOLARE
DIATHERMY UNITS FOR MONOPOLAR AND BIPOLAR SURGERY

MPC System ADC System Output error control

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The Alsatom MPC system is a new generation of ESU (diathermy units) developed using the digital technology with control system by microprocessor, the latest today available for this kind of devices and the only one able to guarantee both performances really suitable to the actual surgical needs and the very high levels of safety, today compulsory according to the latest safety rule for this type of devices (IEC 60601-2-2 III ed., today already European Standard EN) Alsatom line is available with 5 models, 4 for monopolar/bipolar use, 1 for bipolar use only.

Alsatom SU 50-MPC is the basic model for monopolar/bipolar use, has a max output power of 80W, can be used only by single foot-switch and it is suitable for general surgery, ORL dermatology, dental surgery and ect.

Alsatom SU 100-MPC is the same to the model SU 50-MPC but, having a max output power of 140W, it is able to meet better every surgical need, including gynaecology, veterinary and ect. This model is also suitable for fulguration in the bladder under liquid.

Alsatom SU 140-MPC has the same output power of the model SU 100-MPC, controls, complete of specific display, for the independent setting of cutting and coagulating power and it can be used not only by the foot-switch but also by the hand-switch (cutting/coagulating) pencil.

Alsatom SU 140/D-MPC is the most complete model of series. It is the same to the model SU 140-MPC, can be used with handswitch pencil, but it has a technical feature not in common with the units for minor surgery. In fact it is equipped with a double foot-switch for the activation of cutting and coagulation and therefore it results not only the most complete and functional answer for all the above-mentioned needs but even the ideal solution for the first aid, the gastroenterology and ect.

Alsatom SU 140/BD-MPC is the model for bipolar use for cutting or coagulation with a max output power of 100W. It is the same to the model SU 140-MPC with double foot-switch and 2 separate controls, complete of specific display, for the setting of cutting and coagulation powers.

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CONSTRUCTIONAL TECHNOLOGY

● Microprocessor control - MPC system

MPC system is a sophisticated control system of the functioning parameters able to grant a level of performances and safety otherwise unobtainable. The system which performs a complete autodiagnosis both switching on the unit and during the use, allows an extremely safe and easy use because, in case of faults, stops automatically the output power informing the users through specific error codes.

● Automatic adjustment of output powers with dynamic control and stop in case of faults - ADC system with OUTPUT ERROR CONTROL

The output powers are controlled by the microprocessor that ensures an automatic very efficient and dynamic optimisation in every surgical field, carrying out 7000 tests/second thus granting the best surgical results always using the lowest possible power.

This system grants also the maximum safety because, as nowadays it is compulsory according to the latest safety rule, stops in real time the output power, in case of internal faults, and therefore prevents from a sudden and potentially dangerous output power exceeding the selected value.

This system checks also the HF leakage currents and, if there are any load unmatching impedances (patient tissue) reduces automatically the RMS voltage to safety levels.

● Power setting "Micro Type" with control display

Control display for mod. SU140, SU140-D, SU140-BD

● Neutral plate safety circuit and activation signals "sounds and lights".

In compliance with IEC Standard 60601-2-2 III ed. (see technical features)

● Pneumatic foot-switches

All the foot switches are "high safety type" water proof and explosion proof, being pneumatic without electrical current.

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PERFORMANCES

The Alsatom have 5 currents, 3 the mod. 140/BD (bipolar only), to meet every need and all the currents have a starting impulse electronically controlled to get the best results by using the lowest possible output power.

The currents are:

- **CUT:** without modulation for pure cut without haemostasis;
- **BLEND:** modulated combination of currents for cut and haemostasis for cut coagulating
- **COAG:** at high voltage "fulguration" for coagulation without contact and with strong superficial effect. It is indispensable for many surgical needs
- **MICRO:** at low voltage "pin-point, soft", extremely precise with excellent depth effect (ex. for depilation in dermatology);
- **BIPOLAR:** without modulation and extremely precise "micro".

The model Alsatom SU140/BD-MPC has in addition to bipolar coagulation 2 currents for bipolar cut:

- **CUT:** without modulation, for cutting;
- **BLEND:** modulated for cut coagulating

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THE UNITS, THE STANDARD SET OF ACCESSORIES

- B700** - ALSATOM SU 50-MPC, without accessories
B710 - ALSATOM SU 100-MPC, without accessories
B720 - ALSATOM SU 140-MPC, without accessories
B700/A - SET OF STANDARD ACCESSORIES FOR:
ALSATOM SU 50-MPC, SU 100-MPC, SU 140-MPC as follows:
1STOP/PN - single switch pedal, pneumatic, water-proof, explosion-proof
1MPE/F - autoclavable pencil with cable m 2.5
2E1 - straight knife electrode (blade)
1E5 - thick needle electrode
1E7 - thin needle electrode
1E12 - ball electrode Ø 2.5 mm
1E14 - ball electrode Ø 4 mm
1EIP/9 - stainless steel neutral plate with cable 2.5 m
1FFE - fixing belt for neutral plate
B700/B - ditto, with flexible, conductive silicone neutral plate NP/GP
B700/D - ditto, but for dental use, without EIP/9 and SEL/VI, with: EIP/S - hand cylinder neutral plate with cable 2.5 m SEL/D - set of 8 dental electrodes

For manual activation (only mod. SU 140-MPC), see: MPE/CMS, MP/CM

- B730** - ALSATOM SU 140/D-MPC, without accessories
B730/A - SET OF STANDARD ACCESSORIES FOR ALSATOM SU 140/D-MPC identical to B700/A, but with dual switch pedal D-STOP/P
B730/B - ditto, but with flexible, conductive silicone neutral plate NP/GP
B730/D - ditto, but for dental use, without EIP/9 and SEL/VI, with: EIP/S - hand cylinder neutral plate with cable 2.5 m SEL/D - set of 8 dental electrodes

For manual activation see: MPE/CMS, MP/CM (*)

- B740** - ALSATOM SU 140/BD-MPC, with dual switch pedal D-STOP/P

For accessories see: "ACCESSORIES AND CABLES FOR BICOAGULATION" MICROBICOAGULATION, BICOAGULATION OF THE TURBINAL OR LARYNX"

To connect cables of bipolar accessories (electrodes for coagulation or cutting, even for laparoscopic procedures) produced by other manufacturers, see:

"ADAPTERS FOR CONNECTING MONOPOLAR CABLES OR BIPOLAR ACCESSORIES"

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GENERAL DETAIL OF ACCESSORIES

PENCILS

- MPE/F** - autoclavable pencil with cable mt. 2.5
MPE/CMS - hand-switch pencil (CUT/COAG), 100 times autoclavable with cable mt. 3
MP/CM - Ditto, disposable (with knife or needle or ball electrode) (*)

Reusable NEUTRAL ELECTRODES

- EIP/9** - stainless steel neutral electrode with cable mt. 2.5
NP/GP - Ditto, flexible conductive silicone
EIP/S - Ditto, hand cylinder

Disposable NEUTRAL ELECTRODES

- CMS/E** - reusable connecting cable, mt. 2.5
EIP/DA - adhesive electrode for adults (conf. 25 pcs.) (*)
EIP/DP - Ditto, pediatric (*)

TROLLEYS AND FOOT-SWITCHES

- H10** - two-shelves trolley (BxDxH cm 41x36x35)
H10/F - Ditto, with an arm for the tubing of smoke suction unit
STOP/PN - pneumatic, waterproof, explosion proof single foot-switch
D-STOP/P - Ditto, double

SMOKE SUCTION UNITS

- CS200** - smoke suction unit, with filter

PENCILS AND NEEDLES FOR MICRO SURGERY

- MLD/F** - autoclavable light pencil, cable mt. 2
MID - needle adaptor (for pencils MPE/F, MPE/CMS and ect.)
AID - insulated needle - Ø 0.10 mm for depilation
SAD - set of long and non insulated needles Ø 0.10 mm for fine cut
SAD/1 - Ditto, Ø 0.15 mm for fine cut, telangiectasia and ect.
SAD/2 - Ditto, Ø 0.20 mm for fine cut, face depilation, telangiectasia of legs and ect.
SAD/3 - Ditto, Ø 0.40 mm, for fine cut, microcoagulation and ect.

ACCESSORIES AND CABLES FOR BICOAGULATION, MICROBICOAGULATION, BICOAGULATION OF THE TURBINAL OR OF LARYNX, see Alsa catalogue: "ACCESSORIES FOR BICOAGULATION".

ELECTRODES FOR ELETTROCHIRURGIA, FOR GYNECOLOGY AND CONIZATION reusable or disposable, FOR DENTAL USE and ect., see Alsa catalogue: "MONOPOLAR ELECTRODES".

ADAPTERS FOR MONOPOLAR CONNECTING CABLES (ex. Units for endoscopy) OR OF BIPOLAR ACCESSORIES (made by other manufacturers), see Alsa catalogue: "ADAPTERS".

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