

Gynaecological Workplace

reduce both spatial demands in the surgery and the physician's time



Gracie Office

A modern and friendly environment for both the physician and patient emphasizing ergonomics and effectiveness of gynaecological investigation.



All in one place

Gynaecological examination chair Integrated digital video colposcope LCD monitor Sophisticated and complex instrument tray








Video Colposcope An innovative solution provided by BORCAD – integration of a colposcope

in the Gracie chair reduces spatial demands and increase's convenience for the physician during the examination.



Technical parameters

Magnification: 1.5x - 20x optical, 10x digital Visible area diameter: 3 cm (magnification 10x)

Working distance: 24-33 cm

Focus depth: 1 cm (magnification 10x)

LCD diagonal: 17" Resolution: 1280x1024

Light source: LED

Vessel imaging filter: green color

Data formats for storage: png, avi

Output: USB

Photograph resolution: 640*480 px Video resolution: 640*480 px





Technical parameters

Length: 600 mm Height: 1125 mm Width: 520 mm

Waste basket volume: 17 / 13 l Used speculums container volume: 25 l Upper container for speculums volume: 13 l





- GRACIE GYNAECOLOGICAL WORKPLACE
 examination chair with integrated lying position, digital colposcope and LCD
 instrument tray inc. baskets for waste and speculums
- doctor's chair



Designed by Jiří Španihel



Accessories





GKB-054 Shallow instrumental bowl





PVC seat cover



Rotary instrumental bowl (L)



GKB-77 PVC underseat cover



GKB-079 Large bowl



Colposcope holders



GKB-00.01 PVC pedal cover - pair




GKB-021 Doctor's foot support



GKB-044 Socket terminal (DE)



Doctor's mobile chair





GKB-068.X Seat heating

GKB-070 Backup battery



Gynaecological Workplace

