


Cryosurgical Systems: Fast, Precise Freeze and Thaw Control

OUR ONE- AND TWO-TRIGGER SYSTEMS PROVIDE PRECISE CONTROL AND A LARGE SELECTION OF HIGHLY CONDUCTIVE, REUSABLE CRYOTIPS AND DISPOSABLE PLASTIC SHIELDS ARE AVAILABLE



LL100™ System (N₂O) Two-trigger Design and Twice the Freeze Power of Other Systems

Wallach's exclusive two-trigger cryosurgical system delivers fast and precise freeze/defrost control and allows cryotips to be changed at any time without venting.

- More effective freeze/defrost in less time
- Freeze and defrost only when triggered
- Lightweight, well-balanced, and maneuverable
- Use with Wallach's large selection of cryotips and disposable plastic shields

LLCO₂™ System Single-trigger Design is Optimized for CO₂ Gas

The LLCO₂ single-trigger cryosurgical system maximizes temperature and minimizes the clogging and sputtering commonly associated with CO₂ gas.

- Squeeze trigger to freeze; release to defrost
- Use with medical-grade CO₂ gas for optimal performance
- Use with Wallach's standard cryotips with disposable plastic shields

Both the LL100 and the LLCO₂ are available in three-tip package specials for additional savings. Choose three cryotips from our comprehensive offerings when ordering. Wallach recommends the T-1905, T-2514, and T-0500 for GYN procedures and the T-0524, T-0219, and T-1900 for general/family practice procedures. Call Wallach customer service for more information.



ITEM NUMBER	DESCRIPTION
900001	LL100 Cryosurgical System – N ₂ O (Freezer Only) ^{Δ†}
900160	LLCO ₂ Cryosurgical System – CO ₂ (Freezer Only) ^Δ
901060	N ₂ O Cylinder (unfilled)
901061	CO ₂ Cylinder (unfilled)
901011	Cart for 20 lb. Cylinder

WA1000B™ Cryo Console System

Uses LL100™ two-trigger or LLCO₂™ single-trigger freezers

- High-performance features and easy-to-use controls
- Temperature Indicator with integrated timer for precise documentation
- Cart carries a single 20 lb. cylinder
- Use with Wallach's large selection of cryotips and shields

ITEM NUMBER	DESCRIPTION
900506-1-N2O	WA1000B Cryosurgical Console – Blank N ₂ O LL100 ^{Δ†}
900506-1-CO2	WA1000B Cryosurgical Console – Blank CO ₂ LL100 [†]
900506-2-N2O	WA1000B Cryosurgical Console – Digital Temp. (N ₂ O) Indicator, Timer, Base, and LL100 ^Δ
900506-2-CO2	WA1000B Cryosurgical Console – Digital Temp. (CO ₂) Indicator, Timer, Base, and LLCO ₂



Cryosurgical Tips for Use with LL100™ and LLCO₂™

Highly conductive, thermal, reusable, sterilizable cryotips with disposable plastic shields.

 900300AA (T-0200) ^{Δ†}	 900301A (T-0500) ^{Δ†}	 900302A (T-0700) ^{Δ†}	 900303AA (T-0900) ^{Δ†}	 900305AA (T-1300) ^{Δ†}	 900306AA (T-1500) ^{Δ†}
 900200AA (T-0219) ^{Δ†}	 900201AA (T-0219) ^{Δ†}	 900202AA (T-0524) ^{Δ†}	 900203AA (T-0519) ^{Δ†}	 900204AA (T-0832) ^{Δ†}	 900205AA (T-0819) ^{Δ†}
 900206AA (T-0826) ^{Δ†}	 900207AA (T-0823) ^{Δ†}	 900208AA (T-1900) ^{Δ†}	 900209AA (T-1905) ^{Δ†}	 900210AA (T-1920) ^{Δ†}	 900211AA (T-1910) ^{Δ†}
 900212AA (T-1920) ^{Δ†}	 900213AA (T-2500) ^{Δ†}	 900214AA (T-2507) ^{Δ†}	 900215AA (T-2514) ^{Δ†}	 900217AA (T-0832) ^{Δ†}	Disposable Plastic Shields for Cryotips (Bag of 100) 900221AA ^{Δ†}

